

2021 生醫年會優惠方案: 購買 norgen 任一商品(不限下列商品)就送小米手環 5, 滿 1 萬送小米小愛音箱 or 遊戲片 x1(NS/PS4/PS5), 滿 2 萬以上送 Apple HomePod mini, 滿 6 萬以上送 Switch 一台+任一片遊戲片!

贈品只能擇 1 且不能搭配其他優惠方案



Total RNA Extraction and Purification (including microRNA)

Kit Name	Size	Cat. #
Total RNA Purification Kit	50 preps	17200
	100 preps	37500
Total RNA Purification Plus Kit	50 preps	48300
	100 preps	48400
Total RNA Purification 96-Well Kit	2 x 96-well	24300
Total RNA Purification Maxi Kit	8 preps	26800
Total RNA Purification Micro Kit	50 preps	35300
Total RNA Purification Plus Kit	50 preps	48500
microRNA Purification Kit	25 preps	21300
Single Cell RNA Purification Kit	50 preps	51800
Animal Tissue RNA Purification Kit	50 preps	25700
Fatty Tissue RNA Purification Kit	25 preps	36200

Isolate RNA, microRNA, DNA and proteins from the Same Sample

RNA/DNA Purification Kit	50 preps	48700
RNA/DNA Purification Micro Kit	50 preps	50300
RNA/Protein Purification Kit	50 preps	48200
RNA/DNA/Protein Purification Plus Kit	50 preps	47700
RNA/DNA/Protein Purification Plus Micro Kit	50 preps	51600
RNA/DNA/Protein Purification 96-Well Kit	1 x 96-well	51700

Clean-Up and Concentration of microRNA and Total RNA

RNA Clean-Up and Concentration Kit	50 preps	23600
	100 preps	43200
RNA Clean-Up and Concentration 96-Well Kit	2 x 96-well	25100
RNA Clean-Up and Concentration Micro-Elute Kit	50 preps	61000
CleanAll RNA/DNA Clean-Up and Concentration Kit	50 preps	23800
Oligo Clean-Up and Concentration Kit	50 preps	34100
<i>In Vitro</i> Transcription Clean-Up and Concentration Micro Kit	50 preps	52900
<i>In Vitro</i> Transcription Clean-Up and Concentration Mini Kit	50 preps	53000

Stool Purification

Stool Total RNA Purification Kit	50 preps	49500
Stool Nucleic Acid Isolation Kit	50 preps	45600

Cytoplasmic and Nuclear RNA and microRNA Purification

Cytoplasmic and Nuclear RNA Purification Kit	50 preps	21000
	100 preps	27400

FFPE RNA and microRNA Extraction and Purification

FFPE RNA Purification Kit	50 preps	25300
FFPE RNA Purification 96-Well Kit	2 x 96-well	25400
FFPE RNA/DNA Purification Plus Kit	50 preps	54300

Plant and Environmental RNA and microRNA Purification

Kit Name	Size	Cat. #
Plant/Fungi Total RNA Purification Kit	50 preps	25800
	100 preps	31350
Plant/Fungi Total RNA Purification 96-Well Kit	2 x 96-well	31900
Plant microRNA Purification Kit	50 preps	54700
Plant RNA/DNA Purification Kit	50 preps	24400
Soil Total RNA Purification Kit	50 preps	27750
Water RNA/DNA Purification Kit (0.22 µm)	50 preps	26400
Water RNA/DNA Purification Kit (0.45 µm)	50 preps	26450

Urine RNA and microRNA Isolation and Purification

Urine Cell-Free Circulating RNA Purification Mini Kit	50 preps	56900
Urine Cell-Free Circulating RNA Purification Midi Kit	20 preps	57000
Urine Cell-Free Circulating RNA Purification Maxi Kit	10 preps	57100
Urine Cell-Free Circulating and Viral Nucleic Acid Purification Mini Kit	50 preps	56900
	20 preps	57000
Urine Cell-Free Circulating and Viral Nucleic Acid Purification Midi Kit	10 preps	57100
	50 preps	22550
Urine Exfoliated Cell RNA Purification Kit	50 preps	22550
Urine microRNA Purification Kit	50 preps	29000
Urine Exosome RNA Isolation Kit	50 preps	47200

Blood/Plasma/Serum Sample Preparation Kits

Plasma/Serum RNA Purification Mini Kit	50 preps	55000
Plasma/Serum RNA Purification Midi Kit	20 preps	56100
Plasma/Serum RNA Purification Maxi Kit	10 preps	56200
Plasma/Serum RNA/DNA Purification Mini Kit	50 preps	55200
Plasma/Serum Cell-Free Circulating and Viral Nucleic Acid Purification Mini Kit	50 preps	56300
	20 preps	56400
Plasma/Serum Cell-Free Circulating and Viral Nucleic Acid Purification Midi Kit	10 preps	56500
	50 preps	21200
Leukocyte RNA Purification	2 x 96-well	37800
Leukocyte RNA Purification 96-Well Kit	48 preps	52620
Norgen's Blood RNA Preservation and Purification System	96 preps	52720
	48 preps	52600
Preserved Blood RNA Purification Kit (for use with Norgen Blood RNA Preservative Tubes)	96 preps	52700
	50 preps	43400
Preserved Blood RNA Purification Kit I (for use with Tempus Blood RNA Tubes)	50 preps	43400
Preserved Blood RNA Purification Kit II (for use with PAXgene Blood RNA Tubes)	50 preps	43500

浩盟生技有限公司
電話: 02-29083971

只有 norgen RNA spin column 能留住 total RNA 和高純度



左圖為各家廠牌的 RNA kit 中使用 spin column 後, 抽完 RNA 跑膠比較圖

NORGEN
BIOTEK CORP.
The Sample Preparation Experts

HMB
Harmony Biosolution Inc.