

## WizPure™ qPCR Master (SYBR)

產品特色：確保試劑盒提供快速、高特異性和超靈敏的 real-time PCR

- ◇ 即用型解決方案
- ◇ 簡單的反應設置
- ◇ 最低優化要求
- ◇ 抗體介導的熱啟動 PCR
- ◇ 用於通用 PCR 預混液
- ◇ 適用於大多數 real-time 定量 PCR 儀
- ◇ 應用獨特的增強和穩定技術


Step	( ° C)	Time	Cycle
Initial Denaturation	95	5min	1
Denature	95	10-30sec	30-40
Anneal	55-68	10-60sec	
Melting curve analysis	65-95	2-5sec/step	1

Cat No.	Size	Remark
W1711	1ml	100rxn
W1711-5	5x1ml	500rxn
W1711R	1ml	100rxn with ROX dye
W1711R-4	4x1ml	400rxn with ROX dye
W1711RL	1ml	Low ROX dye
W1711RL-5	5x1ml	Low ROX dye
W1711RH	1ml	High ROX dye
W1711RH-5	5x1ml	High ROX dye




➤ 有興趣的同學歡迎索取 FREE SAMPLE

### Amplification



### Melt Peak



wizpure qPCR master (SYBR) 和競爭對手的 real-time PCR 試劑盒在人  $\beta$ -肌動蛋白基因擴增中的比較。

